



PRODUCT SPECIFICATIONS

P/N TOILET SENSOR

1. CONSTRUCTION

- 1.A.** Length of Switch: 24 inches
- 1.B.** Width of Switch: 21 inches
- 1.C.** Top Material: Embossed vinyl, 0.70 inch thick.
- 1.D.** Bottom Material: PVC coated polyester, pebble finish, no-slip
- 1.E.** Cable: 2 x 24 AWG, PVC sheath. Various lengths and connectors are available. Optionally a customer specified wireless transmitter may be installed.

2. OPERATION

- 2.A.** Activation Force: 10 lbs.
- 2.B.** Contact Closure: Normally open, momentary contact switch.
- 2.C.** Max current: 50 mA
- 2.D.** Max voltage. 24 volts.

3. INSTALLATION

- 3.A.** Place toilet sensor in front of toilet, with "ears" on each side of toilet base.

4. CARE

- 4.A.** Cleaning. Wipe down with moistened cloth.
DO NOT MACHINE WASH!

5. DIMENSIONS

